

Quantos Command Set for Quantos System



METTLER TOLEDO

Table of Contents

1	Commands and Responses	3
	QRA 20 – Activate/deactivate user input of value/text.....	3
	QRA 49 – Activate/deactivate message window.....	4
	QRA 61 1 – DO command powder/liquid dosing.....	6
	QRA 61 3 – DO command cut label.....	7
	QRA 61 4 – DO command stop dosing.....	8
	QRD 2 3 7 – Enquiry of the front door position.....	9
	QRD 2 3 8 – Enquiry of the sampler position.....	10
	QRD 2 2 8 – Get sampler status.....	11
	QRD 2 2 9 – Get weighing pan status.....	12
	QRD 2 4 11 – Get head data.....	13
	QRD 2 4 12 – Get sample data.....	15
	QRD 2 5 12 – Get label sample data.....	17
	QRD 2 6 12 – Get protocol sample data.....	18
	QRA 60 2 – Move dosing head pin.....	19
	QRA 60 7 – Move frontdoor.....	20
	QRA 60 8 – Move autosampler position.....	21
	QRD 1 1 1 – Set value tapping before dosing.....	22
	QRD 1 1 2 – Set value tapping while dosing.....	23
	QRD 1 1 3 – Set tapper intensity.....	24
	QRD 1 1 4 – Set tapper duration.....	25
	QRD 1 1 5 – Set target value in mg.....	26
	QRD 1 1 6 – Set target tolerance value in percent.....	27
	QRD 1 1 7 – Set tolerance mode.....	28
	QRD 1 1 8 – Set value sample ID.....	29
	QRD 1 1 9 – Set value pan.....	30
	QRD 1 1 13 – Set value user ID.....	31
	QRD 1 1 14 – Set algorithm selection.....	32
	QRD 1 1 15 – Set AntiStatic Kit.....	33

1 Commands and Responses

QRA_20 – Activate/deactivate user input of value/text

Description

Use QRA_20 to activate or deactivate user input of value or text from the Quantos terminal.

Syntax

Command

QRA_20_x1_ "<text1>"_ "<text2>"_ "<text3>"	Displays an user input dialog.
QRA_20_0	Deactivate active user input of value/text.

First Responses

QRA_20_B	Command executed, user input follows.
QRA_20_A_ "<User_input>"	The user inputted a value that is returned after pressing "OK".
QRA_20_C	QRA_20 command has been aborted by "C" button.
QRA_20_I	Command understood, but not executable at the moment (e.g. there is already an active QRA_20). No second response follows.
QRA_20_L	Command understood, parameter wrong. No second response follows.

Further Responses

QRA_20_A	Command executed, i.e. the last QRA_20 is aborted, i.e. the QRA_20 dialog window is closed.
QRA_20_I	Command understood, but not executable at the moment (e.g. no QRA_20 dialog window is open).

Parameters

Name	Type	Values	Meaning
x1	Integer	8	Alphanumeric input is allowed.
"<text1>"	String		Window title (max. 20 characters).
"<text2>"	String		Text/value to be displayed as default, and to be overwritten by user input (max. 20 characters).
"<text3>"	String		Unit (max. 20 characters).

Comments

- If x1=8, "Text3" can be up to 20 characters. The window title will become "Text1 [Text3]".
- "Home"-, "User Profiles"-, "Select Application"- and "Application Menu"-keys are disabled when softkeys are displayed.

Examples

↓	QRA_20_8_ "Sample_ID"_ "Sample1"_ "-"	Display input dialog to enter the sample ID.
↑	QRA_20_B	Displaying input box.
↑	QRA_20_A_ "Test"	User as entered "Test" as sample ID.

QRA 49 – Activate/deactivate message window

Description

Use QRA_49 to activate or deactivate the message window.

Syntax

Command

QRA_49_x1_x2_ "<text>"	Activates the message window.
QRA_49_0	Closes an open message window.

First Responses

QRA_49_B	Command executed.
QRA_49_A_1	"OK" button is pressed. QRA_49 window will also close.
QRA_49_A_2	"C" button is pressed. QRA_49 window will also close.
QRA_49_I	Command understood, but not executable at the moment.
QRA_49_L	Command understood, parameter wrong, i.e. text is too long or one of the characters cannot be displayed.

Further Responses

QRA_49_A	Command executed, i.e. the last QRA_49 is aborted, i.e. the QRA_49 window is closed.
QRA_49_I	Command understood, but not executable at the moment (e.g. no QRA_49QRA_49 window is open).

Parameters

Name	Type	Values	Meaning
x1	Integer	1	Window without any button.
		2	Window with "OK" button only.
		3	Window with "C" button only.
		4	Window with "OK" and "C" button.
x2	Integer	1	Info icon.
		2	Alert icon.
		3	Stop icon.
		4	Question icon.
		5	X-mark icon.
		6	Date/Time icon.
		7	Sound icon.
		8	Alphabetic icon.
		9	Numeric icon.
"<text>"	String		Sandglass icon.
			Message in the window (max. 240 characters).

Comments

- "Home"-, "User Profiles"-, "Select Application"- and "Application Menu"-keys are disabled when QRA_49 window is displayed. The window has to be closed before they can be used again.
- If no button is visible, the window can be closed using the side door opening buttons.

- If there is an active QRA_49, sending another QRA_49 command will close the current window before displaying a new window with the new setting.

Example

↓	QRA_49_3_1_"This_is_a_sample_text_window."	Activates a message window containing the text "This is a sample text window.".
↑	QRA_49_B	Command executed. This text becomes visible in the window. The window includes "C" button and info icon.
↑	QRA_49_A_2	"C" key is pressed.

QRA 61 1 – DO command powder/liquid dosing

Description

This command starts instantly a powder or liquid dosing.

Syntax

Command

QRA_61_1	Starts dosing.
----------	----------------

Responses

QRA_61_1_B	Executing command.
QRA_61_1_A	Command executed and complete.
QRA_61_1_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRA_61_L	Command understood but not executable (incorrect parameter).

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		6	Weight not stable.
		7	Powderflow error.
		8	Stopped by external action.
		9	Safepos error.
		10	Head not allowed.
		11	Head limit reached.
		12	Head expiry date reached.
		13	Sampler blocked.

Comment

- This dosing uses all the parameters previously set, such as: target, tolerance and powder dosing algorithm.

Example

↓	QRA_61_1	Command to start dosing.
↑	QRA_61_1_B	Executing command.
↑	QRA_61_1_A	Dosing is finished.

QRA 61 3 – DO command cut label

Description

This command cuts a label that has been printed.

Syntax

Command

QRA_61_3	Cuts the label.
----------	-----------------

Responses

QRA_61_3_B	Executing command.
QRA_61_3_A	Command executed and complete.
QRA_61_3_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRA_61_L	Command understood but not executable (incorrect parameter).

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.

Comment

- This command cuts the label if it has been printed through the Quantos command set.

Example

↓	QRA_61	Command to cut label.
↑	QRA_61_3_B	Executing command.
↑	QRA_61_3_A	Command executed and complete.

QRA 61 4 – DO command stop dosing

Description

This command stops the dosing process.

Syntax

Command

QRA_61_4	Stops a running dosing process.
----------	---------------------------------

Responses

QRA_61_4_B	Executing command.
QRA_61_4_A	Command executed and complete.
QRA_61_4_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRA_61_L	Command understood but not executable (incorrect parameter).

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		6	Weight not stable.
		7	Powderflow error.
		8	Stopped by external action.
		9	Safepos error.
		10	Head not allowed.
		11	Head limit reached.
		12	Head expiry date reached.
		13	Sampler blocked.

Example

↓	QRA_61_4	Stops dosing.
↑	QRA_61_4_B	Executing command.
↑	QRA_61_4_A	Dosing process stopped.

QRD 2 3 7 – Enquiry of the front door position

Description

This command returns the position of the Quantos front door.

Syntax

Command

QRD_2_3_7	Requests the front door position.
-----------	-----------------------------------

Responses

QRD_2_3_7_<FrontdoorPosition>_A	Command executed and returns front door position.
QRD_2_3_7_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_3_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<FrontdoorPosition>	Integer	2	Close position.
		3	Open position.
		8	Not detectable.
		9	Running.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- If the front door is unmounted in the system, this command is not executable.

Example

↓	QRD_2_3_7	Requests the front door position.
↑	QRD_2_3_7_3_A	The front door is in the open position.

QRD 2 3 8 – Enquiry of the sampler position

Description

This command returns the current position of the Quantos autosampler QS30.

Syntax

Command

QRD_2_3_8	Requests the sampler position.
-----------	--------------------------------

Responses

QRD_2_3_8_<SamplerPosition>_A	Command executed and returns sampler position.
QRD_2_3_8_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_3_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<SamplerPosition>	Integer	0	Home position.
		1-30	Position.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Returns QRD_2_3_I when sampler is deactivated or not available.

Example

↓	QRD_2_3_8	Requests the sampler position.
↑	QRD_2_3_8_3_A	Sampler is on position 3.

QRD 2 2 8 – Get sampler status

Description

Requests the activation status of the Quantos autosampler QS30.

Syntax

Command

QRD_2_2_8	Requests the sampler status.
-----------	------------------------------

Responses

QRD_2_2_8_<SamplerStatus>_A	Command executed and complete.
QRD_2_2_8_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_2_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<SamplerStatus>	Integer	0	Sampler is switched off/disabled.
		1	Sampler is switched on/selected.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Check sampler status before using it to avoid errors.

Example

↓	QRD_2_2_8	Requests the sampler status.
↑	QRD_2_2_8_1_A	Sampler is mounted and activated.

QRD 2 2 9 – Get weighing pan status

Description

This command detects whether the weighing pan is empty or contains a vial. Initialise the weighing pan before using this function by calling the set empty pan function, **see** QRD_1_1_9.

Syntax

Command

QRD_2_2_9	Requests the weighing pan status.
-----------	-----------------------------------

Responses

QRD_2_2_9<PanStatus>_A	Command executed and complete.
QRD_2_2_9_I<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_2_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<PanStatus>	Integer	0	Pan empty.
		1	Pan not empty.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Be sure to initialise the empty weighing pan before using this command, **see** QRD_1_1_9.

Example

↓	QRD_2_2_9	Requests the weighing pan status.
↑	QRD_2_2_9_1_A	Weighing pan is not empty (there is a vial placed).

QRD 2 4 11 – Get head data

Description

This command reads the content of the dose head RFID chip and returns it as an XML stream.

Syntax

Command

QRD_2_4_11	Get dose head data.
------------	---------------------

Responses

QRD_2_4_11_B	Command understood, sending xml.
<XML_Dataset>	xml result stream.
QRD_2_4_11_A	Command executed and complete.
QRD_2_4_11_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_4_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Tag	Type	Values	Meaning
<XML_Dataset>	<?xml version="1.0" encoding= "ISO-8 859-1"?> <Info_head xmlns:xsi= "http:// www.w3.org/ 2001/XMLSchema instance" xsi:no NamespaceSchema Location=" Info_head.xsd"> <Timestamp>	String	Text	Time stamp.
	<Substance>	String	Text	Substance name.
	<Lot_ID>	String	Text	Lot number of the substance.
	<User_ID>	String	Text	User ID of the currently logged in user.
	<Filling_date>	String	Text	Filling date of the dose head.
	<Exp._date>	String	Text	Expiry date of the substance.
	<Retest_date>	String	Text	Retest date of the substance.
	<Content Unit="mg">	String	Text	Fill weight of powder.
	<Rem._dosages>	String	Text	Remaining number of doses.
	<Terminal_SNR>	String	Text	Serial number of the terminal.
	<Bridge_SNR>	String	Text	Serial number of the weighing bridge.
	<Balance_Type>	String	Text	Balance type.
	<Balance_ID>	String	Text	Balance ID.
	<Last_cal.>	String	Text	Date of last calibration.
	<Option_SNR>	String	Text	Serial number of the ethernet option.
	<Dose_unit_SNR>	String	Text	Serial number of the dosing unit.
	<Appl._Name>	String	Text	Fixed application name.
	<Var1>	String	Text	User variable 1.

	<Var2>	String	Text	User variable 2.
	<Var3>	String	Text	User variable 3.
	<Var4>	String	Text	User variable 4.
	<Label>	String	Text	Description of a user variable.
	<Value>	String	Text	Content of the user variable.
	<Title_1>	String	Text	Title 1.
	<Title_2>	String	Text	Title 2.
	<Date_Time>	String	Text	Date and time stamp.
	<Levelcontrol>	String	Text	Level control setting.
	<Head_prod._ date>	String	Text	Production date of dose head.
	<Head_type>	String	Text	Dose head type.
	<Head_ID>	String	Text	Unique ID of the dose head.
	<Dose_limit>	String	Text	Dose head limit.
	<Accuracy Unit="%">	String	Text	Accuracy (not used).
	<Dosing_ counter>	String	Text	Number of doses done with this head.
	<Rem._quantity Unit="mg">	String	Text	Remaining amount of powder in the dose head.
	</Info_head>	String	Text	-

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Powder and liquid dose heads have slightly different output files, see the xml output specification.

Example

↓	QRD_2_4_11	Get dose head xml data.
↑	QRD_2_4_11_B	Sending xml data.
↑	QRD_2_4_11_A	XML sent completely.

QRD 2 4 12 – Get sample data

Description

Receives the result data from the last dispense (powder or liquid).

Syntax

Command

QRD_2_4_12	Get dose sample data.
------------	-----------------------

Responses

QRD_2_4_12_B	Command understood, sending xml.
<XML_Dataset>	xml result stream.
QRD_2_4_12_A	Command executed and complete.
QRD_2_4_12_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_4_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Tag	Type	Values	Meaning
<XML_Dataset>	<?xml version="1.0" encoding= "ISO-8 859-1"?> <Info_head xmlns:xsi= "http:// www.w3.org/ 2001/XMLSchema instance" xsi:no NamespaceSchema Location=" Info_head.xsd"> <Timestamp>	String	Text	Time stamp.
	<Substance>	String	Text	Substance name.
	<Lot_ID>	String	Text	Lot number of the substance.
	<User_ID>	String	Text	User ID of the currently logged in user.
	<Filling_date>	String	Text	Filling date of the dose head.
	<Exp._date>	String	Text	Expiry date of the substance.
	<Retest_date>	String	Text	Retest date of the substance.
	<Content Unit="mg">	String	Text	Fill weight of powder.
	<Rem._dosages>	String	Text	Remaining number of doses.
	<Terminal_SNR>	String	Text	Serial number of the terminal.
	<Bridge_SNR>	String	Text	Serial number of the weighing bridge.
	<Balance_Type>	String	Text	Balance type.
	<Balance_ID>	String	Text	Balance ID.
	<Last_cal.>	String	Text	Date of last calibration.
	<Option_SNR>	String	Text	Serial number of the ethernet option.
	<Dose_unit_SNR>	String	Text	Serial number of the dosing unit.
	<Appl._Name>	String	Text	Fixed application name.
	<Var1>	String	Text	User variable 1.

	<Var2>	String	Text	User variable 2.
	<Var3>	String	Text	User variable 3.
	<Var4>	String	Text	User variable 4.
	<Label>	String	Text	Description of a user variable.
	<Value>	String	Text	Content of the user variable.
	<Title_1>	String	Text	Title 1.
	<Title_2>	String	Text	Title 2.
	<Date_Time>	String	Text	Date and time stamp.
	<Levelcontrol>	String	Text	Level control setting.
	<Head_prod._ date>	String	Text	Production date of dose head.
	<Head_type>	String	Text	Dose head type.
	<Head_ID>	String	Text	Unique ID of the dose head.
	<Dose_limit>	String	Text	Dose head limit.
	<Accuracy Unit="%">	String	Text	Accuracy (not used).
	<Dosing_ counter>	String	Text	Number of doses done with this head.
	<Rem._quantity Unit="mg">	String	Text	Remaining amount of powder in the dose head.
	</Info_head>	String	Text	-

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- The output files from liquid and solid dispenses may differ.

Example

↓	QRD_2_4_12	Get sample xml data.
↑	QRD_2_4_12_B	Sending xml data.
↑	QRD_2_4_12_A	XML sent completely.

QRD 2 5 12 – Get label sample data

Description

This command triggers the Quantos label printing function for a sample label.

Syntax

Command

QRD_2_5_12	Triggers printing of a sample label.
------------	--------------------------------------

Responses

QRD_2_5_12_B	Command understood, sending data to the label printer.
QRD_2_5_12_<Labelprinter>_A	Command executed and complete.
QRD_2_5_12_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_5_L	Command understood but not executable (incorrect parameter).

Parameter

Name	Type	Values	Meaning
<Labelprinter>	String		Print out the sample data to the label printer.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- This command works both with powder and liquid dispenses. The output may differ according to the configuration in the Quantos application settings.

Example

↓	QRD_2_5_12	Send sample data to the label printer.
↑	QRD_2_5_12_B	Printing started.
↑	QRD_2_5_12_A	Printing completely sent.

QRD 2 6 12 – Get protocol sample data

Description

This command triggers the printing of the sample label through Quantos using the strip printer.

Syntax

Command

QRD_2_6_12	Prints sample protocol.
------------	-------------------------

Responses

QRD_2_6_12_B	Command understood, sending protocol data to the strip printer.
QRD_2_6_12_A	Command executed and complete.
QRD_2_6_12_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_2_6_L	Command understood but not executable (incorrect parameter).

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- This command works both with powder and liquid dispenses. The output may differ according to the configuration in the Quantos application settings.

Example

↓	QRD_2_6_12	Trigger printing a sample protocol.
↑	QRD_2_6_12_B	Printing started.
↑	QRD_2_6_12_A	Printing completely sent.

QRA 60 2 – Move dosing head pin

Description

This command unlocks or locks the dosing head to the dosing unit.

Syntax

Command

QRA_60_2_OpeningPosition>	Locks or unlocks dosing head.
---------------------------	-------------------------------

Responses

QRA_60_2_B	Executing command.
QRA_60_2_A	Command executed and complete.
QRA_60_2_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRA_60_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<OpeningPosition>	Integer	3	Unlock position.
		4	Lock position.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.

Example

↓	QRA_60_2_3	Command to unlock.
↑	QRA_60_2_B	Executing command.
↑	QRA_60_2_A	Dosing head unlocked.

QRA 60 7 – Move frontdoor

Description

Opens or closes the Quantos front door to allow easy access to the dose head.

Syntax

Command

QRA_60_7<FrontdoorPosition>	Open or close the front door.
-----------------------------	-------------------------------

Responses

QRA_60_7_B	Executing command.
QRA_60_7_A	Command executed and complete.
QRA_60_7_I_<Error_code>	Command understood but currently not executable, according to.
QRA_60_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<FrontdoorPosition>	Integer	2	Close position.
		3	Open position.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Also check locking and unlocking the dose head, **see** QRA_60_2.

Example

↓	QRA_60_7_3	Open the front door.
↑	QRA_60_7_B	Quantos is opening the front door (busy).
↑	QRA_60_7_A	Action completed, front door is open.

QRA 60 8 – Move autosampler position

Description

This command moves the autosampler to the selected position.

Syntax

Command

QRA_60_8_<AutosamplerPosition>	Set the autosampler position.
--------------------------------	-------------------------------

Responses

QRA_60_8_B	Executing command.
QRA_60_8_A	Command executed and complete.
QRA_60_8_I_<Error_code>	Command understood but currently not executable, according to.
QRA_60_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<AutosamplerPosition>	Integer	0	Home position.
		1 ... 30	Respective position.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		8	Stopped by external action.
		13	Sampler blocked.

Comment

- The autosampler configuration in the terminal must be considered to this command. The configuration set the smallest moving step of the autosampler.

Example

↓	QRA_60_8_3	Moves autosampler to position 3.
↑	QRA_60_8_B	Executing command.
↑	QRA_60_8_A	Command executed and complete.

QRD 1 1 1 – Set value tapping before dosing

Description

This command activates or deactivates the tapping before dosing function.

Syntax

Command

QRD_1_1_1_<TapBeforeDosing>	Sets the tap before dosing function.
-----------------------------	--------------------------------------

Responses

QRD_1_1_1_A	Setting stored.
QRD_1_1_1_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<TapBeforeDosing>	Integer	0	Tap before dosing deactivated.
		1	Tap before dosing activated.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Tapping before dosing can help to dispense powders. This command is parameterised by the tapper intensity and duration.

Example

↓	QRD_1_1_1_1	Set tap before dosing to ON.
↑	QRD_1_1_1_A	Setting stored.

QRD 1 1 2 – Set value tapping while dosing

Description

This command activates or deactivates the tap while dosing function.

Syntax

Command

QRD_1_1_2_<TapWhileDosing>	Activates or deactivates the tap while dosing function.
----------------------------	---

Responses

QRD_1_1_2_A	Setting stored.
QRD_1_1_2_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<TapWhileDosing>	Integer	0	Deactivate tapping while dosing.
		1	Activate tapping while dosing.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- The actual usage of tap while dosing is controlled from the algorithm. Tapper intensity or duration has no impact on this. Note that tapping can lead to compacting effects with some powders.

Example

↓	QRD_1_1_2_1	Activate tap while dosing function.
↑	QRD_1_1_2_A	Setting stored.

QRD 1 1 3 – Set taper intensity

Description

Sets the intensity of the taper.

Syntax

Command

QRD_1_1_3_<Intensity>	Set the taper intensity.
-----------------------	--------------------------

Responses

QRD_1_1_3_A	Intensity setting stored.
QRD_1_1_3_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_3_L	Command understood but not executable (incorrect parameter).

Parameter

Name	Type	Values	Meaning
<Intensity>	Integer	10 - 100	Intensity of the tap before dosing function in percent.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- The intensity value only applies on tap before dosing or manual tapping.

Example

↓	QRD_1_1_3_100	Set taper intensity to 100%.
↑	QRD_1_1_3_A	Setting stored.

QRD 1 1 4 – Set taper duration

Description

This command specifies the duration in seconds of the tap before dosing.

Syntax

Command

QRD_1_1_4_<Duration>	Set the taper duration.
----------------------	-------------------------

Responses

QRD_1_1_4_A	Setting stored.
QRD_1_1_4_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameter

Name	Type	Values	Meaning
<Duration>	Integer	1 - 10	Duration in seconds.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- The duration value is only applied for tap before dosing or manual tapping.

Example

↓	QRD_1_1_4_2	Set tap duration to 2 seconds.
↑	QRD_1_1_4_A	Setting stored.

QRD 1 1 5 – Set target value in mg

Description

Sets the target amount for liquid and solid doses in milligrams.

Syntax

Command

QRD_1_1_5_<Target>	Set the target value in mg.
--------------------	-----------------------------

Responses

QRD_1_1_5_A	Setting stored.
QRD_1_1_5_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameter

Name	Type	Values	Meaning
<Target>	Float	0.10 – 250000.00	Enter the target value of substance to be dosed in mg with two decimal places.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- Note that both for liquid and solid, the target is specified in milligrams. This value should be set before executing a dose, **see** QRA_61_1.

Example

↓	QRD_1_1_5_50.00	Set target value to 50.00 mg.
↑	QRD_1_1_A	Setting stored.

QRD 1 1 6 – Set target tolerance value in percent

Description

Defines the tolerance value for powder dispensing in percent.

Syntax

Command

QRD_1_1_6_<Tolerance>	Set the target tolerance value in percent.
-----------------------	--

Responses

QRD_1_1_6_A	Setting stored
QRD_1_1_6_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameter

Name	Type	Values	Meaning
<Tolerance>	Float	0.1 – 40.0	Enter the target tolerance value of substance to be dosed in percent with one decimal place.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comments

- The tolerance setting only affects the powder dispensing. Note that the tolerance mode, **see** QRD_1_1_7 affects the way the tolerance is handled.
- It is recommended to have a tolerance of (absolute) 0.5 mg in order to achieve repeatable dispensing results.

Example

↓	QRD_1_1_6_1.0	Set the target tolerance value to 1.0%.
↑	QRD_1_1_6_A	Setting stored.

QRD 1 1 7 – Set tolerance mode

Description

This command defines the operation of the tolerance. 0/+ tolerance mode only allows overdosing, +/- tolerance mode applies the tolerance band below and above the target value.

Syntax

Command

QRD_1_1_7_<ToleranceMode>	Set value of the tolerance mode.
---------------------------	----------------------------------

Responses

QRD_1_1_7_A	Setting stored.
QRD_1_1_7_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<ToleranceMode>	Integer	0	+/- tolerance mode.
		1	0/+ tolerance mode.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	2	Another job is running.
		3	Timeout.
		5	Not allowed at the moment.
		8	Stopped by external action.

Comment

- The tolerance mode is used in conjunction with the tolerance setting, **see** QRD_1_1_6 for powder dispensing.

Example

↓	QRD_1_1_7_0	Set tolerance mode to +/- tolerance.
↑	QRD_1_1_7_A	Tolerance mode set.

QRD 1 1 8 – Set value sample ID

Description

This command sets an ID for a sample.

Syntax

Command

QRD_1_1_8_<SampleID>	Sets value of sample ID.
----------------------	--------------------------

Responses

QRD_1_1_8_A	Command executed and complete.
QRD_1_1_8_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<SampleID>	String	Max. 20 chars	Enter the sample ID.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.

Comment

- The sample ID set remains the default as long the sample ID question is excluded from the dosing steps.

Example

↓	QRD_1_1_8_ID1	Set sample ID with value "ID1".
↑	QRD_1_1_8_A	Sample ID set.

QRD 1 1 9 – Set value pan

Description

This command defines whether the pan is empty.

Syntax

Command

QRD_1_1_9_<Pan>	Set weighing pan as empty.
-----------------	----------------------------

Responses

QRD_1_1_9_B	Executing command.
QRD_1_1_9_A	Command executed and complete.
QRD_1_1_9_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<Pan>	Integer	0	Pan is empty.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.
		6	Weight not stable.

Example

↓	QRD_1_1_9_0	Set weighing pan as empty.
↑	QRD_1_1_9_A	Empty weighing pan set.

QRD 1 1 13 – Set value user ID

Description

This command sets an ID for the user.

Syntax

Command

QRD_1_1_13_<UserID>	Set value of user ID.
---------------------	-----------------------

Responses

QRD_1_1_13_A	Command executed and complete.
QRD_1_1_13_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<UserID>	String	Max. 20 chars	Enter the user ID.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.

Comment

- The user ID set becomes the default which can be changed during the dosing steps.

Example

↓	QRD_1_1_13_User1	Set user ID with name "User1".
↑	QRD_1_1_13_A	User ID set.

QRD 1 1 14 – Set algorithm selection

Description

This command selects the powder dosing algorithm.

Syntax

Command

QRD_1_1_14_<AlgoSelection>	Set value of the algorithm selection.
----------------------------	---------------------------------------

Responses

QRD_1_1_14_A	Setting stored.
QRD_1_1_14_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_I	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<AlgoSelection>	Integer	0	Select algorithm "Standard".
		1	Select algorithm "Advanced".

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.

Comment

- These algorithms only influence the powder dispense.

Example

↓	QRD_1_1_14_0	Set algorithm.
↑	QRD_1_1_14_A	Algorithm selection set.

QRD 1 1 15 – Set AntiStatic Kit

Description

This command activates or deactivates the AntiStatic Kit.

Syntax

Command

QRD_1_1_15_<AntiStaticKit>	Sets the AntiStatic Kit.
----------------------------	--------------------------

Responses

QRD_1_1_15_A	Setting stored.
QRD_1_1_15_I_<Error_code>	Command understood but currently not executable, according to the following error code.
QRD_1_1_L	Command understood but not executable (incorrect parameter).

Parameters

Name	Type	Values	Meaning
<AntiStaticKit>	Integer	0	AntiStatic Kit deactivated.
		1	AntiStatic Kit activated.

Error Codes

Name	Type	Values	Meaning
<Error_code>	Integer	1	Not mounted.
		2	Another job is running.
		3	Timeout.
		4	Not selected.
		5	Not allowed at the moment.

Comment

- This command does not recognize whether the AntiStatic Kit is physically connected or not.

Example

↓	QRD_1_1_15_1	Set AntiStatic Kit ON.
↑	QRD_1_1_15_A	Setting stored.



Good Weighing Practice™

GWP® is the global weighing standard, ensuring consistent accuracy of weighing processes, applicable to all equipment from any manufacturer. It helps to:

- Choose the appropriate balance or scale
- Calibrate and operate your weighing equipment with security
- Comply with quality and compliance standards in laboratory and manufacturing

► www.mt.com/GWP

www.mt.com/quantos

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contact

Subject to technical changes.
© Mettler-Toledo GmbH 08/2016
30348483A en



30348483